


Some classes of Dangerous Chemicals





Pyrophorics

a substance that will ignite spontaneously in air or upon exposure to water or water vapor

Organolithiums: t-BuLi, n-BuLi, MeLi

Organozincs: Et₂Zn

Organomagnesiums: EtMgBr

Organoaluminums: Me₃Al

Metal powders: Zn, Al, Mg, Ti, Ni, Pd

Alkali metals: K, Na **Metal hydrides:** LiAlH₄

Non-metal alkyls: BR₃, PR₃, AsR₃




Often used in conjunction with flammables!

Working Safely with Pyrophorics

- 1.) Don't work alone!
 - notify others, remind them of safety procedures
 - beware of the location of the class D fire extinguisher
- 2.) Wear the right protection
 - nitrile gloves or (nitrile + neoprene) gloves
 - goggles or face shield
 - buttoned-up fire-resistant lab coat
- 3.) Check the state of the reagent and your work space
 - corroded septum?
 - excess flammable material or clutter in hood?
- 4.) Use the proper technique
 - syringe or cannula?

Syringe and Cannulation Techniques

Watch UCLA pyrophoric safety video: http://www.youtube.com/watch?v=RaMXwNBAbxc



Carcinogens

- 1.) Label known carcinogens as such.
- 2.) Store in carcinogen cabinet or in secondary containment.
- 3.) Change gloves after using carcinogens.

Some common carcinogens:

Explosives: Azides

What NOT to do with NaN₃:

- use with chlorinated solvents (di- and tri- azidomethanes)
- mix with CS2, bromine, Bronsted acids, or heavy metals.

In general, "safe" azides conform to the formula:

(#O + #C)/(#N) > or = 3

What to avoid with organic azides:

- concentrating to dryness
- putting under vacuum
- ground-glass joints
- mechanical shock
- distillation
- exposure to light and silica

What to definitely do:

- put in separate waste container
- avoid mixing with acid

Peroxides and peroxide formers

What to avoid:

- metals
- heat, friction, shock
- vacuum

How to hydrolize:

- Use 80 parts water, 20 parts NaOH, and small amount of surfactant
- Slowly pour peroxide in solution with stirring (use ~10X more alkaline solution)

How to test for peroxides:

Use starch-iodide paper.

- immerse in solution for 1 sec
- breathe on strip for 30 sec

No peroxide



Common peroxide formers:

Corrosives: Aqua Regia and Piranha

AVOID CONTACT WITH ORGANIC SOLVENTS

Aqua regia:

- 1:3 nitric to hydrochloric acid
- use proper ventilation!
 HNO₃ + 3 HCl → NOCl + Cl₂ + 2 H₂O
- don't use is closed containers
- dilute with water slowly and neutralize

Piranha:

- 1:3 30% hydrogen peroxide to sulfuric acid
- mix slowly!
- dilute with water and neutralize

Try other cleaning methods first.

Other Common Dangerous Chemicals

Reproductive toxins

- carbon disulfide
- toluene
- pyridine
- HF

Acute toxins

- diazomethane
- sodium azide
- OsO₄
- organomercury
- organotin